



UMS
UNIVERSITI MALAYSIA SABAH



FACULTY OF SCIENCE AND NATURAL RESOURCES

SV40303

SCIENTIFIC DATA VISUALISATION

SEMESTER 1, 2025/ 2026 SESSION

PREPARED FOR:

PROF. ABDULLAH BADE

ASSIGNMENT 2: Advanced Medical Data Visualizer (DICOM)

PREPARED BY:

NO.	NAME	NO. MATRICS
1.	WARDAH HANAANI BINTI AKHMAL	BS22110420
2.	NUR AINA NADIRAH BINTI KAMARUL AZHAR	BS22110050

TABLE OF CONTENT

No.	Content	Page
1.	System Introduction	3
2.	System Architecture	4
3.	System Flowchart	5
4.	Algorithms	6 – 12
5.	Output	13 – 17

Introduction

For Assignment 2, we were tasked to prepare an advanced medical visualizer that handles DICOM (Digital Imaging and Communications in Medical) dataset. Several classes are introduced with functions that help in handling the dataset to be visualize.

DetachableSliceViewer class is the class used to handle the window that shows the slicing of the DICOM dataset loaded into the system. There are a few functions used to show the slicing of the DICOM dataset loaded, such as the functions that update the color changed on any tab functions that are rendered on the 3D image, orientation, transformation, and the slider that shows which part of the slicing is currently shown on the window.


HelpDialog class is used to show a dialog using the QDialog that has the information and direction on how to use or where to click to make changes to the 3D image that will or has been loaded into the system.

Class myVTK is used to render the VTK-related variables, such as the renderer, window, interactor, mapper, actors, and scene, and functions that read the DICOM dataset, set the coordinate from the DICOM loaded with the coordinate, change the transformations such as translation, rotation, and scaling, and other rendering pipelines. There are also functions for camera settings, the slicing camera that shows different orientations of the slicing and information on the slicing to be rendered on the slicing window.


MeasurementInteractor class is used to handle all measurement-related functions. There are several functions that handle the calculations for distance, angle, area, and volume such as `_measure_distance`, `_measure_angle`, `_measure_area`, and `_measure_vollume` with a few additional functions that help to finalize the points marked for area and volume as they did not have an assigned number of points needed to measure. To visualize the calculations, functions `_mark_point`, `_draw_line`, `_draw_polygon`, and `_draw_tetrahedron` are used to show the line and points for marking distance between two points and angle between three points, and lines, points, and polygon or tetrahedron drawn to visualize the area and volume. A function named `on_left_click` is used to control the marking of the points on the 3D image shown using mouse.

The AdvancedDICOMViewer class is where the user interface (UI) is designed by setting the UI design style, creating the left panel and tabs with different functions. Then, there are several functions used to connect the toggles in the tabs and functions that connect the input from the toggles by appending the functions from other classes to make changes to the 3D image.

System Architecture



System Flowchart



Algorithms

Algorithm 1 Main Application Controller

```
1: CLASS AdvancedDICOMViewer EXTENDS QMainWindow
2:                                     ▷ MVC Controller coordinating all components
3:
4: ATTRIBUTES:
5:   vtk_app: myVTK                                     ▷ 3D rendering engine reference
6:   loaded_datasets: Dict[str, DatasetInfo]             ▷ Loaded DICOM datasets
7:   slice_viewer: DetachableSliceViewer                ▷ 2D slice viewer dock
8:   tabs: QTabWidget                                   ▷ 7 functional control tabs
9:   vtk_widget: QVTKRenderWindowInteractor             ▷ 3D viewport
10:
11: procedure SETUP_UI
12:   Create central widget and main layout
13:   Create splitter(LeftPanel, RightPanel)
14:   LeftPanel ← create_left_panel()                   ▷ Control tabs
15:   RightPanel ← create_viewport_panel()               ▷ 3D view
16:   Setup menus, toolbar, statusbar
17: end procedure
18:
19: procedure INITIALIZE_VTK
20:   vtk_app ← new myVTK()
21:   vtk_app.start(vtk_widget)                         ▷ Initialize VTK pipeline
22: end procedure
```

Algorithm 2 VTK Rendering Engine (myVTK)

```
1: CLASS myVTK
2:                                     ▷ Core rendering engine managing VTK pipeline
3:
4: ATTRIBUTES:
5:   renderer: vtkRenderer                 ▷ 3D scene renderer
6:   window: vtkRenderWindow               ▷ Render window
7:   interactor: vtkRenderWindowInteractor ▷ User interaction
8:   imageData: vtkImageData              ▷ DICOM volume data
9:   volumeMapper: vtkSmartVolumeMapper    ▷ Volume rendering
10:  volumeProperty: vtkVolumeProperty     ▷ Visual properties
11:
12: procedure LOAD_DICOM_DATA(directory)
13:   reader ← new vtkDICOMImageReader()
14:   reader.SetDirectoryName(directory)
15:   reader.Update()
16:   imageData ← reader.GetOutput()
17:                                     ▷ Extract metadata: dimensions, spacing, scalar range
18:   dims ← imageData.GetDimensions()
19:   spacing ← imageData.GetSpacing()
20:   scalar_range ← imageData.GetScalarRange()
21:   return success
22: end procedure
23:
24: procedure SETUP_RENDERING_PIPELINE(vtk_widget)
25:   renderer ← new vtkRenderer()
26:   window ← vtk_widget.GetRenderWindow()
27:   window.AddRenderer(renderer)
28:   interactor ← window.GetInteractor()
29:   interactor.SetInteractorStyle(vtkInteractorStyleTrackballCamera)
30: end procedure
```

Algorithm 3 DICOM Data Processing Pipeline

Require: DICOM directory containing slice images

Ensure: 3D rendered volume with interactive controls

```
1:
2: procedure DICOM PROCESSING PIPELINE(directory)
3:   Phase 1: Data Loading & Reconstruction
4:     Read DICOM slices → vtkDICOMImageReader
5:     Stack slices → vtkImageData (3D volume)
6:     Extract metadata → dimensions, spacing, HU range
7:
8:   Phase 2: Transformation Pipeline
9:     transformation ← new vtkTransform()
10:    reslice ← new vtkImageReslice()
11:    reslice.SetInputData(imageData)
12:    reslice.SetResliceTransform(transformation)
13:    transformed_data ← reslice.GetOutput()
14:
15:   Phase 3: Visualization Paths
16:     Path A: 3D Volume Rendering
17:       volumeMapper.SetInputData(transformed_data)
18:       Setup transfer functions (color, opacity)
19:       volume ← new vtkVolume()
20:       renderer.AddVolume(volume)
21:
22:     Path B: 2D Slice Extraction
23:       sliceMapper ← new vtkImageSliceMapper()
24:       sliceMapper.SetInputData(transformed_data)
25:       sliceMapper.SetOrientation(axial/coronal/sagittal)
26:       imageSlice ← new vtkImageSlice()
27:       renderer.AddViewProp(imageSlice)
28:
29:   Phase 4: Interactive Rendering
30:     renderer.ResetCamera()
31:     window.Render()
32:     interactor.Initialize()
33: end procedure
```

Algorithm 4 Interactive Measurement System

```
1: CLASS MeasurementInteractor EXTENDS vtkInteractorStyleTrackballCamera
2:                                     ▷ Handles interactive measurements in 3D space
3:
4: ATTRIBUTES:
5:   points: List[(x,y,z)]                                     ▷ Measurement points
6:   mode: {Distance, Angle, Area, Volume}
7:   drawn_actors: List[vtkActor]                               ▷ Visual feedback
8:
9: procedure ON_LEFT_CLICK(obj, event)
10:   click_pos ← GetInteractor().GetEventPosition()
11:   picker ← new vtkCellPicker()
12:   picker.SetTolerance(0.001)
13:   picker.Pick(click_pos[0], click_pos[1], 0, renderer)
14:   picked_pos ← picker.GetPickPosition()
15:
16:   if picked_pos ≠ (0,0,0) then
17:     points.append(picked_pos)
18:     _mark_point(picked_pos)                                   ▷ Add visual marker
19:     _process_measurement_mode()
20:   end if
21: end procedure
22:
23: procedure _MEASURE_DISTANCE
Require: len(points) = 2
24:   p1, p2 ← points[0], points[1]
25:   dx, dy, dz ← p2 - p1
26:   distance ←  $\sqrt{dx^2 + dy^2 + dz^2}$ 
27:   spacing ← imageData.GetSpacing()
28:   distance_mm ← distance × mean(spacing)
29:   _draw_line(p1, p2)                                           ▷ Visual feedback
30:   _update_ui(f"Distance: {distance_mm:.2f} mm")
31:   points.clear()
32: end procedure
33:
34: procedure _MEASURE_ANGLE
Require: len(points) = 3
35:   p1, p2, p3 ← points[0], points[1], points[2]
36:   v1 ← [p1[i] - p2[i] for i in 0..2]
37:   v2 ← [p3[i] - p2[i] for i in 0..2]
38:   dot ←  $\sum_{i=0}^2 v1[i] \times v2[i]$ 
39:   mag1 ←  $\sqrt{\sum v1[i]^2}$ 
40:   mag2 ←  $\sqrt{\sum v2[i]^2}$ 
41:   angle ←  $\arccos(\text{dot} / (\text{mag1} \times \text{mag2}))$ 
42:   angle_deg ← angle × 180/π
43:   _update_ui(f"Angle: {angle_deg:.2f}°")
44:   points.clear()
45: end procedure
```

Algorithm 5 Organ Segmentation System

Require: vtkImageData (DICOM volume), organ_type**Ensure:** 3D surface model of segmented organ

```
1:
2: procedure SEGMENT_ORGAN(min_hu, max_hu, organ_name)
3:   Step 1: Thresholding
4:   threshold  $\leftarrow$  new vtkImageThreshold()
5:   threshold.SetInputData(imageData)
6:   threshold.ThresholdBetween(min_hu, max_hu)
7:   threshold.SetInValue(1)  $\triangleright$  Foreground
8:   threshold.SetOutValue(0)  $\triangleright$  Background
9:   threshold.Update()
10:
11:   Step 2: Connectivity Filtering
12:   connectivity  $\leftarrow$  new vtkImageConnectivityFilter()
13:   connectivity.SetInputConnection(threshold.GetOutputPort())
14:   connectivity.SetExtractionModeToLargestRegion()
15:   connectivity.Update()
16:   binary_mask  $\leftarrow$  connectivity.GetOutput()
17:
18:   Step 3: Surface Generation
19:   marching  $\leftarrow$  new vtkMarchingCubes()
20:   marching.SetInputData(binary_mask)
21:   marching.SetValue(0, 0.5)  $\triangleright$  Iso-surface value
22:   marching.ComputeNormalsOn()
23:   marching.Update()
24:   surface  $\leftarrow$  marching.GetOutput()
25:
26:   Step 4: Visualization
27:   mapper  $\leftarrow$  new vtkPolyDataMapper()
28:   mapper.SetInputConnection(marching.GetOutputPort())
29:   mapper.ScalarVisibilityOff()
30:
31:   actor  $\leftarrow$  new vtkActor()
32:   actor.SetMapper(mapper)
33:   actor.GetProperty().SetColor(0.9, 0.7, 0.3)  $\triangleright$  Gold color
34:   actor.GetProperty().SetOpacity(0.9)
35:
36:   renderer.AddActor(actor)
37:   volume.SetVisibility(False)  $\triangleright$  Hide main volume
38:   window.Render()
39: end procedure
40:  $\triangleright$  Organ HU Ranges (Hounsfield Units)
41: Bone: [200, 3000]
42: Lungs: [-1000, -600]
43: Liver: [50, 70]
44: Kidneys: [30, 50]
45: Soft Tissue: [-50, 100]
46: Muscle: [10, 40]
```

Algorithm 6 Inter-module Communication

```
1: Communication Pattern: Event-Driven with Signals/Slots
2:
3: UI ↔ VTK Engine:
4: UI.Event → Controller.Method → VTK.Update → Render
5: VTK.RenderComplete → Callback → UI.StatusUpdate
6:
7: Slice Viewer ↔ Main View Synchronization:
8: procedure SYNC_SLICES
9:   if slice_viewer.sync_enabled then
10:     slice_viewer.sliceChanged → main_view.set_slice
11:     main_view.orientationChanged → slice_viewer.set_orientation
12:   end if
13: end procedure
14:
15: Measurement System Flow:
16: UserClick → Picker → 3D Coordinates
17: → Geometry Calculation → Visual Feedback
18: → UI Display → History Storage
19:
20: Data Export Flow:
21: ExportRequest → Data Processor → Format Converter
22: → File Writer → Status Notification
```

Algorithm 7 Rendering Optimization Strategies

```
1: Optimization Level 1: Adaptive Sampling
2: procedure SET_RENDERING_QUALITY(level) level 1
3:   volumeMapper.SetSampleDistance(4.0) 2
4:   volumeMapper.SetSampleDistance(2.0) 3
5:   volumeMapper.SetSampleDistance(1.0) 4
6:   volumeMapper.SetSampleDistance(0.5) 5
7:   volumeMapper.SetSampleDistance(0.25)
8: end procedure
9:
10: Optimization Level 2: Surface Simplification
11: procedure OPTIMIZE_SURFACE(quality_level)
12:   decimate ← new vtkDecimatePro()
13:   reduction ← [0.9, 0.7, 0.4, 0.2, 0.0][quality_level-1]
14:   decimate.SetTargetReduction(reduction)
15:   decimate.PreserveTopologyOn()
16: end procedure
17:
18: Memory Management:
19: • Lazy loading of DICOM slices
20: • vtkSmartPointer for automatic memory management
21: • Cache frequently used transformations
22: • Progressive rendering for large datasets
23:
24: UI Responsiveness:
25: • Asynchronous data loading
26: • Debounced UI updates for sliders
27: • Separate rendering thread
```

```

1: STRUCT ApplicationState
2:                                     ▷ Encapsulates complete application state
3:
4: Datasets:
5:   active_dataset: DatasetInfo
6:   loaded_datasets: Map[id: DatasetInfo]
7:
8: Visualization Settings:
9:   view_mode: {3D, AxialSlice, CoronalSlice, SagittalSlice}
10:  render_mode: {Volume, Surface, MIP, MinIP}
11:  transfer_function: Preset or Custom
12:  window/level: (width, center)
13:
14: Transformations:
15:   translation: (x, y, z) in mm
16:   rotation: (rx, ry, rz) in degrees
17:   scaling: factor
18:
19: Measurements:
20:   active_measurements: List[Measurement]
21:   history: List[MeasurementResult]
22:
23: Segmentation:
24:   segmentations: Map[name: SegmentationMask]
25:   active_segmentation: SegmentationMask
26:
27: Camera State:
28:   position: (x, y, z)
29:   focal_point: (x, y, z)
30:   view_up: (x, y, z)
31:   clipping_range: (near, far)

```

Algorithm 8 Application State Representation

Algorithm 9 Multi-format Export Pipeline

```

1: PIPELINE ExportData(format, data_type)
2: format DICOM
3: writer ← new vtkDICOMWriter()
4: writer.SetInputData(imageData)
5: writer.SetDirectoryName(output_dir)
6: writer.SetFilePrefix("modified.")
7: writer.Write() PNG/JPEG/TIFF
8: window_to_image ← new vtkWindowToImageFilter()
9: window_to_image.SetInput(window)
10: writer ← format_specific_writer()
11: writer.SetFileName(output_file)
12: writer.Write() VTK Volume
13: writer ← new vtkXMLImageDataWriter() ▷ .vti
14: writer.SetInputData(imageData)
15: writer.SetFileName(output_file)
16: writer.Write() Complete Scene
17: multiBlock ← new vtkMultiBlockDataSet()
18: multiBlock.SetBlock(0, imageData)
19: multiBlock.SetBlock(1, settings.json)
20: writer ← new vtkXMLMultiBlockDataWriter()
21: writer.Write() ▷ .vtm file
22: =0

```

Pattern	Implementation	Purpose
MVC	AdvancedDICOMViewer (Controller)	Separation of concerns
Observer	PyQt5 Signals/Slots	Event-driven updates
Strategy	Rendering modes, Measurement algorithms	Algorithm interchangeability
Composite	Multi-block scene organization	Hierarchical data management
Factory	Dataset loaders, Export writers	Object creation abstraction
Singleton	Application instance	Single point of access

Algorithm 10 Application Startup Sequence


```

1: procedure MAIN
2:   app ← new QApplication(sys.argv)
3:   app.setStyle('Fusion')
4:
5:   window ← new AdvancedDICOMViewer()
6:   window.show()
7:
8:                                     ▷ Setup global shortcuts
9:   fullscreen_shortcut ← new QShortcut("F11", window)
10:  fullscreen_shortcut.activated.connect(window.toggle_fullscreen)
11:
12:  escape_shortcut ← new QShortcut("Escape", window)
13:  escape_shortcut.activated.connect(exit_fullscreen_handler)
14:
15:  sys.exit(app.exec_())
16: end procedure


```

Output


Data Tab




View Tab



Color Tab




Segment Tab




Measure Tab


Measuring Distance




Measuring Angle




Measuring Area



Measuring Volume




Transform Tab




Render Tab


Volume Rendering



Surface Rendering



Maximum Intensity Projection



Minimum Intensity Projection

